

Product Datasheet

Anti-Human CD25 Capture Antibody, Azide Free Clone B-G3 879.500.001

Product Name	Anti-Human CD25 Capture Antibody, Azide Free Clone B-G3
Catalog Number	879.500.001
Description	This Antibody can be used as Capture in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	PHA activated T cells
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,5-2,5µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	PHA activated T cells
Specificity	Recognises the Interleukin 2 Receptor a chain (IL-2 Ra), a 55 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells