

Product Datasheet

Anti-Human CD130 Capture Antibody, Azide Free Clone B-P4 879.520.001

Product Name	Anti-Human CD130 Capture Antibody, Azide Free Clone B-P4
Catalog Number	879.520.001
Description	This Antibody can be used as Capture in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Natural soluble gp130
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,5-5µg/ml
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Gp130, common subunit for IL-6, IL-11, OSM, LIF, CNTF, CT-1 receptors, a 130-140 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells