

Product Datasheet

Anti-Human Perforin Detection Antibody, Biotin Conjugated Clone B-D48 879.860.002

Product Name Anti-Human Perforin Detection Antibody, Biotin Conjugated Clone
B-D48

Catalog Number 879.860.002

Description This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay.

Species/Host Mouse

Conjugation Biotin

Tested Applications ELISA, ELISpot

Immunogen Natural human perforin

Brand Diaclone

Form/Appearance Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide

Storage Stable at +2-8°C for 12 months.

Note The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA/ELISpot: 0,05-0,5µg/ml

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises human perforin, a 70kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells