

Product Datasheet

Anti-h ST2 10202 SPTN-5 100681

Product Name	Anti-h ST2 10202 SPTN-5
Catalog Number	100681
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	CLIA, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Recognizes both free sST2 and sST2 bound to IL-33
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Recombinant ST2 antigen, Medix Biochemica, 710020
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.3-6.8
Cross Reactivity	N/D
Specificity	Antibody recognizes human soluble ST2
Shelf Life	Unspecified