

Product Datasheet

Anti-Human IP-10 Detection Antibody, Biotin Conjugated Clone B-C55 879.950.002

Product Name	Anti-Human IP-10 Detection Antibody, Biotin Conjugated Clone B-C55
Catalog Number	879.950.002
Description	This biotinylated Antibody can be used as Detection in a human sandwich Immunoassay.
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant IP-10
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide
Storage	Stable at +2-8°C for 12 months.
Note	The suggested coating concentration range below should be optimised by each laboratory for each application. ELISA: 0,05-0.5µg/ml
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant IP-10
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells