

# Product Datasheet

## Human IFN gamma ELISA Kit 950.000.192

Product Name	Human IFN gamma ELISA Kit
Catalog Number	950.000.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	12.5-400 pg/ml
Sensitivity	5 pg/ml
Tested Applications	ELISA
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g
Specificity	Recognizes both natural and recombinant human IFN-g
Sample Types	Serum, Plasma, Cell culture supernatant