

Product Datasheet

Human IL-2 ELISA Kit 950.010.096

Product Name	Human IL-2 ELISA Kit
Catalog Number	950.010.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	7 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-2.
Specificity	Recognizes both natural and recombinant human IL-2
Sample Types	Serum, Plasma, Cell culture supernatant