

# Product Datasheet

## Human IL-6 ELISA Kit 950.030

Product Name	Human IL-6 ELISA Kit
Catalog Number	950.030
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	6.25-200 pg/ml
Sensitivity	2 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. No interference with soluble receptor gp80.
Specificity	Recognizes both natural and recombinant human IL-6
Sample Types	Serum, Plasma, Cell culture supernatant