

Product Datasheet

Human IL-10 ELISA Kit 950.060

Product Name	Human IL-10 ELISA Kit
Catalog Number	950.060
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	12.5-400 pg/ml
Sensitivity	4.9 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10.
Specificity	Recognizes both natural and recombinant human IL-10
Sample Types	Serum, Plasma, Cell culture supernatant