

## Product Datasheet

### Human IL-10 ELISA Kit 950.060.192

<b>Product Name</b>	Human IL-10 ELISA Kit
<b>Catalog Number</b>	950.060.192
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	12.5-400 pg/ml
<b>Sensitivity</b>	4.9 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	NULL
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10.
<b>Specificity</b>	Recognizes both natural and recombinant human IL-10
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant