

Product Datasheet

Human IL-12p40/p70 ELISA Kit 950.075

Product Name	Human IL-12p40/p70 ELISA Kit
Catalog Number	950.075
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	62.5-2000 pg/ml
Sensitivity	20 pg/ml
Tested Applications	ELISA
Alternative Names	IL-12/IL23p40
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reacts with feline and canine IL-12 p40.
Specificity	Recognizes the p40 monomer, dimer or heterodimer (p40 + p35) forms of the molecule
Sample Types	Serum, Plasma, Cell culture supernatant