

## Product Datasheet

### Anti-h Cystatin C 10002 SPTN-5 100690

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Anti-h Cystatin C 10002 SPTN-5   |
| <b>Catalog Number</b>      | 100690   |
| <b>Description</b>         | Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| <b>Tested Applications</b> | FIA  |
| <b>Brand</b>               | Medix Biochemica   |
| <b>Form/Appearance</b>     | Liquid, may turn slightly opaque during storage  |
| <b>Concentration</b>       | 5.0 mg/ml (+/- 10 %)   |
| <b>Storage</b>             | +2-8°C   |
| <b>Isotype</b>             | IgG2b  |
| <b>Clonality</b>           | Monoclonal   |
| <b>Epitope</b>             | N/D  |
| <b>Purity</b>              | ≥ 95 %   |
| <b>Affinity constant</b>   | N/D  |
| <b>Associated Products</b> | Recombinant Cystatin C antigen, Medix Biochemica, 610100   |
| <b>Buffer</b>              | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative                   |
| <b>IEF Profile</b>         | 6.2-7.3  |
| <b>Cross Reactivity</b>    | Does not recognize recombinant Cystatin D, F, S, SA, or SN.  |
| <b>Specificity</b>         | Antibody recognizes human Cystatin C   |
| <b>Shelf Life</b>          | Unspecified  |