

Product Datasheet

Anti-h Cystatin C 10002 SPTN-5 100690

Product Name	Anti-h Cystatin C 10002 SPTN-5
Catalog Number	100690
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG2b
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Recombinant Cystatin C antigen, Medix Biochemica, 610100
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.2-7.3
Cross Reactivity	Does not recognize recombinant Cystatin D, F, S, SA, or SN.
Specificity	Antibody recognizes human Cystatin C
Shelf Life	Unspecified