

Product Datasheet

Anti-h Inhibin beta A 9504 SPTN-5 100709

Product Name	Anti-h Inhibin beta A 9504 SPTN-5
Catalog Number	100709
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Amino acids 81-112 (CVPTKLRPMSMLYYDDGQNIKKDIQNMIVVEE)
Purity	≥ 95 %
Affinity constant	N/D
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	7.1-8.5
Cross Reactivity	N/D
Specificity	Antibody recognizes human inhibin beta A subunit, also known as activin beta A
Shelf Life	36 months