

Product Datasheet

Anti-Rotavirus VP6 11101 SPTN-5 100715

Product Name	Anti-Rotavirus VP6 11101 SPTN-5
Catalog Number	100715
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Recombinant Rotavirus VP6 antigen, Medix Biochemica 610145
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.8-7.5
Cross Reactivity	N/D
Specificity	Antibody recognizes Rotavirus VP6
Matched Pairs	Forms pair with itself
Shelf Life	18 months