

Product Datasheet

Anti-h KIM-1 10101 SPTN-5 100738

Product Name	Anti-h KIM-1 10101 SPTN-5
Catalog Number	100738
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.2-6.8
Cross Reactivity	Recognizes TIM-3/HAVcr-2/KIM-3. Does not recognize TIM-4/TIMD-4.
Specificity	Antibody recognizes human kidney injury molecule 1
Shelf Life	Unspecified