

Product Datasheet

Anti-h KIM-1 10103 SPTN-5 100740

Product Name	Anti-h KIM-1 10103 SPTN-5
Catalog Number	100740
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.3-7.2
Cross Reactivity	Does not recognize TIM-3/HAVcr-2/KIM-3 or TIM-4/TIMD-4.
Specificity	Antibody recognizes human kidney injury molecule 1
Shelf Life	Unspecified