

Product Datasheet

Anti-h Lp-PLA2 11202 SPTN-5 100751

Product Name	Anti-h Lp-PLA2 11202 SPTN-5
Catalog Number	100751
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Note	Casein is recommended to be used in blocking buffer instead of BSA.
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Associated Products	Recombinant Lipoprotein-associated Phospholipase A2 (Lp-PLA2) antigen 400-60
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.5-7.6
Cross Reactivity	N/D
Specificity	Antibody recognizes human lipoprotein-associated phospholipase A2
Shelf Life	24 months