

## Product Datasheet

### Anti-h AMH 11302 SPTN-5 100757

<b>Product Name</b>	Anti-h AMH 11302 SPTN-5
<b>Catalog Number</b>	100757
<b>Description</b>	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
<b>Tested Applications</b>	LF, CLIA, FIA
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	5.0 mg/ml (+/- 10 %)
<b>Storage</b>	+2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Epitope</b>	N-terminal region of AMH within amino acids Arg26-Arg451
<b>Purity</b>	≥ 95 %
<b>Affinity constant</b>	N/D
<b>Associated Products</b>	<b>Recombinant Anti-Mullerian Hormone (AMH 19-449) Antigen LA578</b>
<b>Buffer</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	5.3-6.0
<b>Cross Reactivity</b>	N/D
<b>Specificity</b>	Antibody recognizes human anti-Müllerian hormone (AMH)
<b>Matched Pairs</b>	FIA: 11301 (capture) - 11303 (detection), 11303 - 11301 and 11309 - 11302 LF: 11304 (membrane) - 11303 (particles) CLIA: 11303 (capture) - 11309 (detection), 11304 - 11309 and 11303 - 11302
<b>Shelf Life</b>	36 months