

Product Datasheet

Anti-h AMH 11304 SPTN-5 100759

Product Name	Anti-h AMH 11304 SPTN-5
Catalog Number	100759
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, CLIA, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N-terminal region of AMH within amino acids Arg26-Arg451
Purity	≥ 95 %
Affinity constant	N/D
Associated Products	Recombinant Anti-Mullerian Hormone (AMH 19-449) Antigen LA578
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	5.8-6.7
Cross Reactivity	N/D
Specificity	Antibody recognizes human anti-Müllerian hormone (AMH)
Matched Pairs	FIA: 11301 (capture) - 11303 (detection), 11303 - 11301 and 11309 - 11302 LF: 11304 (membrane) - 11303 (particles) CLIA: 11303 (capture) - 11309 (detection), 11304 - 11309 and 11303 - 11302
Shelf Life	36 months