

Product Datasheet

Anti-h S100B 11402 SPTN-5 100779

Product Name	Anti-h S100B 11402 SPTN-5
Catalog Number	100779
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Associated Products	Recombinant S100B antigen, Medix Biochemica 710014, 710044 and 710045
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	7.0-7.5
Cross Reactivity	Does not recognize S100A8 or S100A9
Specificity	Antibody recognizes human S100B
Shelf Life	Unspecified