

Product Datasheet

Anti-h D-Dimer 1409 SPTN-5 100800

Product Name	Anti-h D-Dimer 1409 SPTN-5
Catalog Number	100800
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	IT, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	Not Applicable (N/A)
Associated Products	Native D-Dimer, Lee Biosolutions, Cat. 200-09, 200-12 and 200-13
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.1-6.8
Cross Reactivity	Does not recognize human fibrinogen
Specificity	Antibody recognizes human D-dimer
Shelf Life	24 months