

## Product Datasheet

### Anti-Thyroxine 6902 SPTN-5 100801

<b>Product Name</b>	Anti-Thyroxine 6902 SPTN-5
<b>Catalog Number</b>	100801
<b>Description</b>	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
<b>Tested Applications</b>	FIA
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	5.0 mg/ml (+/- 10 %)
<b>Storage</b>	+2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Epitope</b>	N/D
<b>Purity</b>	≥ 95 %
<b>Affinity constant</b>	N/D
<b>Associated Products</b>	Thyroxine (T4) antigen, BSA conjugate, Lee Biosolutions 581-10
<b>Buffer</b>	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	6.7-7.6
<b>Cross Reactivity</b>	T3 (3,3',5-triiodothyronine) 1.5 %, rT3 (3,3',5'-triiodothyronine, reverse T3) 1.5 %, 3,5-diiiodothyronine < 0.1 %
<b>Specificity</b>	Antibody recognizes human thyroxine
<b>Shelf Life</b>	18 months