

Product Datasheet

Anti-h SAA 2212 SPTN-5 100809

Product Name	Anti-h SAA 2212 SPTN-5
Catalog Number	100809
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, IT, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	aa 28-35
Purity	≥ 95 %
Associated Products	Native human SAA, EastCoast Bio cat LA449, recombinant human SAA, Medix Biochemica cat 610070 and Lee Biosolutions cat 527-18
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	7.0 - 7.9
Cross Reactivity	Recognizes human SAA1, human SAA2, feline and equine SAA. Does not recognize human SAA4, canine or bovine SAA.
Specificity	Antibody recognizes human serum amyloid A (SAA)
Shelf Life	24 months