

Product Datasheet

Anti-h Prealbumin 11605 SPTN-5 100832

Product Name	Anti-h Prealbumin 11605 SPTN-5
Catalog Number	100832
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	IT, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Associated Products	Native Prealbumin/Transthyretin antigen 496-10 and 496-11.
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.3-6.7
Cross Reactivity	N/D
Specificity	Antibody recognizes human prealbumin
Shelf Life	Unspecified