

Product Datasheet

Propoxyphene (PPX) Antibody 100895

Product Name	Propoxyphene (PPX) Antibody
Catalog Number	100895
Description	Protein A purified monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, ELISA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	≥ 1.5 mg/ml
Storage	+2-8°C
Note	This product was previously available with Catalog codes P001M and CB302. It has been produced using the same original mouse monoclonal hybridoma cell line.
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	N/D
Buffer	0.05M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN ₃ as a preservative
IEF Profile	6.7 - 7.6
Cross Reactivity	N/D
Specificity	Antibody recognizes Propoxyphene
Matched Pairs	170032 (BSA conjugate as capture)
Shelf Life	24 months