

## Product Datasheet

### 6-Monoacetylmorphine (6-MAM)-BSA Antigen 170031

<b>Product Name</b>	6-Monoacetylmorphine (6-MAM)-BSA Antigen
<b>Catalog Number</b>	170031
<b>Description</b>	Synthetic 6-Monoacetylmorphine, BSA conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Lot-specific
<b>Storage</b>	Short term storage ( $\leq 2$ weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	This product was previously available with Catalog codes P83-99-41A and LA263. It has been produced using the same conjugation protocol.
<b>Buffer</b>	Distilled water
<b>Matched Pairs</b>	140041, 100983, JP091, JP092, JP093, HM445, HM446 or HM461 (detection)
<b>Shelf Life</b>	60 months