

Product Datasheet

Morphine-BSA Antigen 170036

Product Name	Morphine-BSA Antigen
Catalog Number	170036
Description	Synthetic Morphine, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	Lot-specific
Storage	Short term storage (≤ 2 weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P81-99-01A. It has been produced using the same conjugation protocol.
Associated Products	100897 (detection)
Buffer	Distilled water
Shelf Life	60 months