

Product Datasheet

Recombinant Human Interleukin-1 Receptor Accessory Protein IL1RAP (CHO) 705-H23-100

Product Name	Recombinant Human Interleukin-1 Receptor Accessory Protein IL1RAP (CHO)
Catalog Number	705-H23-100
Description	Recombinant Protein
Species/Host	CHO Cells
Alternative Names	IL-1 RAcP (Interleukin-1 Receptor Accessory Protein) , IL1 R3 , C3orf13
Brand	Diaclone
Form/Appearance	Liquid ; Phosphate-buffered saline pH7.4. Sterile-filtered through 0.22 µm.
Storage	Store at minimum -20°C. Avoid repeated freeze-thaw cycles.
Purity	>95% by SDS PAGE
Biological Activity	Recognized by anti-IL-1RAP (B-L43) in ELISA
Sequence	Ser21-Glu359 (ECD)
Tag	HIS Tag C-Terminus