

## Product Datasheet

### Fentanyl-BSA Antigen 170040

<b>Product Name</b>	Fentanyl-BSA Antigen
<b>Catalog Number</b>	170040
<b>Description</b>	Synthetic Fentanyl, BSA conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Lot-specific
<b>Storage</b>	Short term storage ( $\leq 2$ weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	This product was previously available with a Catalog code P81-99-51A. It has been produced using the same conjugation protocol.
<b>Buffer</b>	Distilled water
<b>Matched Pairs</b>	100998, 100974, P01-99-53R-IF, P01-99-54R-IF, P01-99-55R-IF (detection)
<b>Shelf Life</b>	60 months