

## Product Datasheet

### Recombinant Human SARS-CoV-2 Nucleoprotein (HEK) 715-H17-0BU

<b>Product Name</b>	Recombinant Human SARS-CoV-2 Nucleoprotein (HEK)
<b>Catalog Number</b>	715-H17-0BU
<b>Description</b>	Recombinant Protein
<b>Species/Host</b>	HEK cells
<b>Alternative Names</b>	SARS-CoV nucleoprotein , Nucleocapsid protein , N protein
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Liquid PBS
<b>Storage</b>	Store at minimum -20°C. Avoid repeated freeze-thaw cycles
<b>Purity</b>	>95% by SDS page
<b>Biological Activity</b>	Recognized by Covid-19 positive patients serum samples
<b>Sequence</b>	Met1-Ala419
<b>Tag</b>	His Tag C-Terminus