

## Product Datasheet

### Zolpidem-BSA Antigen 170041

<b>Product Name</b>	Zolpidem-BSA Antigen
<b>Catalog Number</b>	170041
<b>Description</b>	Synthetic Zolpidem, BSA conjugated
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	Lot-specific
<b>Storage</b>	Short term storage ( $\leq 2$ weeks): 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	This product was previously available with a Catalog code P83-99-11A. It has been produced using the same conjugation protocol.
<b>Buffer</b>	Distilled water
<b>Matched Pairs</b>	JP046, JP047, JP048, HM311 or HM312 (detection)
<b>Shelf Life</b>	60 months