

Product Datasheet

Recombinant Human SARS-CoV-2 Nucleoprotein (E. Coli) 725-H17-0BU

Product Name	Recombinant Human SARS-CoV-2 Nucleoprotein (E. Coli)
Catalog Number	725-H17-0BU
Description	Recombinant Protein
Species/Host	E. coli
Alternative Names	SARS-CoV nucleoprotein , Nucleocapsid protein , N protein
Brand	Diaclone
Form/Appearance	Liquid Tris HCl 20 mM pH 8.0 ; NaCl 150 mM
Storage	Store at minimum -20°C. Avoid repeated freeze-thaw cycles
Purity	>95% by SDS page
Biological Activity	Recognized by Covid-19 positive patients serum samples
Sequence	Met1-Ala419
Tag	His Tag C-Terminus