

Product Datasheet

Recombinant Pepsinogen I, 100 µg 610000

Product Name	Recombinant Pepsinogen I, 100 µg
Catalog Number	610000
Description	Recombinant human pepsinogen I protein.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	Storage at 2-8 °C. After reconstitution stable for 2 months at 2-8 °C and 2 months at -20 °C.
Associated Products	Anti-h Pepsinogen I 8003, 8009, 8015, and 8016
Buffer	Lyophilized from 50 mM Tris-HCl, pH 7.5; 150 mM NaCl, 0.5 µg/ml pepstatin A; containing 6 % sucrose as a stabilizer
Molecular Weight	Predicted molecular weight: 40 kDa
Shelf Life	36 months