

## Product Datasheet

### Recombinant Pepsinogen I, 100 µg 610000

<b>Product Name</b>	Recombinant Pepsinogen I, 100 µg
<b>Catalog Number</b>	610000
<b>Description</b>	Recombinant human pepsinogen I protein.
<b>Species/Host</b>	E. coli
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Storage at 2-8 °C. After reconstitution stable for 2 months at 2-8 °C and 2 months at -20 °C.
<b>Associated Products</b>	Anti-h Pepsinogen I 8003, 8009, 8015, and 8016
<b>Buffer</b>	Lyophilized from 50 mM Tris-HCl, pH 7.5; 150 mM NaCl, 0.5 µg/ml pepstatin A; containing 6 % sucrose as a stabilizer
<b>Molecular Weight</b>	Predicted molecular weight: 40 kDa
<b>Shelf Life</b>	36 months