

## Product Datasheet

### Human IL-6 ELISpot Kit 856.021.001PC

<b>Product Name</b>	Human IL-6 ELISpot Kit
<b>Catalog Number</b>	856.021.001PC
<b>Description</b>	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	NULL
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. No interference with soluble receptor gp80.
<b>Specificity</b>	Recognizes natural human IL-6