

## Product Datasheet

### Human IL-6 ELISpot Set 856.021.001S

<b>Product Name</b>	Human IL-6 ELISpot Set
<b>Catalog Number</b>	856.021.001S
<b>Description</b>	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 sterile plate.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	NULL
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Note</b>	Spots read visually or with reader
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. No interference with soluble receptor gp80.
<b>Specificity</b>	Recognizes natural human IL-6