

Product Datasheet

Recombinant NGAL antigen, 100 µg 610012

Product Name	Recombinant NGAL antigen, 100 µg
Catalog Number	610012
Description	Recombinant human neutrophil gelatinase-associated lipocalin (NGAL) protein with a C-terminal histidine tag.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Storage	2-8 °C
Associated Products	Anti-h NGAL 4202 Anti-h NGAL 4203 Anti-h NGAL 4204 Anti-h NGAL 4205
Buffer	Lyophilized from phosphate buffered saline, pH 7.4; containing 3 % sucrose, 2 % D-mannitol and 0.01 % Tween 20 as stabilizers
Molecular Weight	Predicted molecular weight: 22 kDa
Shelf Life	Unspecified