
Product Datasheet

Human IL-6 ELISpot Set 856.021.005P

Product Name	Human IL-6 ELISpot Set
Catalog Number	856.021.005P
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 5 non-sterile plates.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human cytokines. No interference with soluble receptor gp80.
Specificity	Recognizes natural human IL-6