

Product Datasheet

Recombinant Myoglobin, 100 µg 610030

| | |
|----------------------------|---|
| Product Name | Recombinant Myoglobin, 100 µg |
| Catalog Number | 610030 |
| Description | Recombinant human myoglobin protein with a C-terminal histidine tag. |
| Species/Host | E. coli |
| Brand | Medix Biochemica |
| Form/Appearance | Lyophilized |
| Storage | Storage at 2-8°C. After reconstitution stable for 1 week at 2-8°C and 1 month at -20°C. |
| Associated Products | Anti-h Myoglobin 7001 Anti-h Myoglobin 7004 Anti-h Myoglobin 7005 |
| Buffer | Lyophilized from 50 mM Tris-HCl, pH 8.0; 150 mM NaCl; containing 3 % sucrose, 2 % D-mannitol and 0.01 % Tween 20 as stabilizers |
| Molecular Weight | Predicted molecular weight: 18 kDa |
| Shelf Life | Unspecified |