

## Product Datasheet

### Recombinant Myoglobin, 100 µg 610030

<b>Product Name</b>	Recombinant Myoglobin, 100 µg
<b>Catalog Number</b>	610030
<b>Description</b>	Recombinant human myoglobin protein with a C-terminal histidine tag.
<b>Species/Host</b>	E. coli
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Storage at 2-8°C. After reconstitution stable for 1 week at 2-8°C and 1 month at -20°C.
<b>Associated Products</b>	Anti-h Myoglobin 7001 Anti-h Myoglobin 7004 Anti-h Myoglobin 7005
<b>Buffer</b>	Lyophilized from 50 mM Tris-HCl, pH 8.0; 150 mM NaCl; containing 3 % sucrose, 2 % D-mannitol and 0.01 % Tween 20 as stabilizers
<b>Molecular Weight</b>	Predicted molecular weight: 18 kDa
<b>Shelf Life</b>	Unspecified