

Product Datasheet

Human IL-10 ELISpot Kit 856.031.001PC

Product Name	Human IL-10 ELISpot Kit
Catalog Number	856.031.001PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10.
Specificity	Recognizes natural human IL-10