

Product Datasheet

Human TNF alpha ELISpot Kit 856.041.001PC

Product Name	Human TNF alpha ELISpot Kit
Catalog Number	856.041.001PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Tested Applications	ELISpot
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with macac TNF-a. No interference with TNFR I and TNFR II
Specificity	Recognizes natural human TNF-a