

## Product Datasheet

### Human IFN gamma ELISpot Kit 856.051.001PC

<b>Product Name</b>	Human IFN gamma ELISpot Kit
<b>Catalog Number</b>	856.051.001PC
<b>Description</b>	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	IFN-g, IFN-gamma
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g.
<b>Specificity</b>	Recognizes natural human IFN-g