

Product Datasheet

Human IFN gamma ELISpot Set 856.051.020S

Product Name	Human IFN gamma ELISpot Set
Catalog Number	856.051.020S
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 20 sterile plates.
Tested Applications	ELISpot
Alternative Names	IFN-g, IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g.
Specificity	Recognizes natural human IFN-g