
Product Datasheet

Human IL-12p70 ELISpot Set 856.081.001P

Product Name	Human IL-12p70 ELISpot Set
Catalog Number	856.081.001P
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 non-sterile plate.
Tested Applications	ELISpot
Alternative Names	IL-12
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes natural human IL-12. Specific to p70 heterodimer.