

## Product Datasheet

### Recombinant Calprotectin antigen, 100 µg 610061

<b>Product Name</b>	Recombinant Calprotectin antigen, 100 µg
<b>Catalog Number</b>	610061
<b>Description</b>	Recombinant human calprotectin heterotetramer (S100A8/S100A9) protein with an N-terminal histidine tag in the S100A8 subunit.
<b>Species/Host</b>	E. coli
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized from 0.2 µm filtered solution
<b>Storage</b>	Storage at 2-8 °C. After reconstitution stable for 1 month at 2-8 °C and 1 month at -20 °C. Reconstituted product is recommended to be stored in plastic (PP) containers.
<b>Associated Products</b>	Anti-h Calprotectin 3403 Anti-h Calprotectin 3404 Anti-h Calprotectin 3405 Anti-h Calprotectin 3406 Anti-h Calprotectin 3407
<b>Buffer</b>	Lyophilized from 20 mM Tris-HCl, 100 mM NaCl, 1 mM EDTA, pH 8.0; containing 6 % sucrose as a stabilizer
<b>Molecular Weight</b>	Predicted molecular weights: S100A8 12 kDa, S100A9 13 kDa
<b>Shelf Life</b>	60 months