

Product Datasheet

Recombinant Calprotectin antigen, 100 µg 610061

Product Name	Recombinant Calprotectin antigen, 100 µg
Catalog Number	610061
Description	Recombinant human calprotectin heterotetramer (S100A8/S100A9) protein with an N-terminal histidine tag in the S100A8 subunit.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized from 0.2 µm filtered solution
Storage	Storage at 2-8 °C. After reconstitution stable for 1 month at 2-8 °C and 1 month at -20 °C. Reconstituted product is recommended to be stored in plastic (PP) containers.
Associated Products	Anti-h Calprotectin 3403 Anti-h Calprotectin 3404 Anti-h Calprotectin 3405 Anti-h Calprotectin 3406 Anti-h Calprotectin 3407
Buffer	Lyophilized from 20 mM Tris-HCl, 100 mM NaCl, 1 mM EDTA, pH 8.0; containing 6 % sucrose as a stabilizer
Molecular Weight	Predicted molecular weights: S100A8 12 kDa, S100A9 13 kDa
Shelf Life	24 months