

Product Datasheet

Human CD178 ELISpot Set 856.091.001P

Product Name	Human CD178 ELISpot Set
Catalog Number	856.091.001P
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 non-sterile plate.
Tested Applications	ELISpot
Alternative Names	FasL, CD95L
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes natural human Fas ligand