

Product Datasheet

Human CD178 ELISpot Kit 856.091.005PC

Product Name	Human CD178 ELISpot Kit
Catalog Number	856.091.005PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Tested Applications	ELISpot
Alternative Names	FasL, CD95L
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes natural human Fas ligand